



[Knowledgebase](#) > [Moto360](#) > [Product Features](#) > [Dexcom App on Wear OS: What You Need to Know](#)

Dexcom App on Wear OS: What You Need to Know

Inas B. - 2025-10-23 - [Product Features](#)

The Dexcom app brings Continuous Glucose Monitoring (CGM) to your Wear OS smartwatch, letting you keep an eye on your glucose levels with ease. To ensure a smooth experience, it's important to understand the app's requirements and limitations.

Wear OS Support:

- Dexcom confirms that the current Wear OS version supports their app, allowing you to view CGM data conveniently on your wrist.

Transmitter Pairing:

- Note that your smartwatch **cannot pair directly with the Dexcom transmitter**. The app works in tandem with a compatible mobile device that communicates with the transmitter.

Mobile Device Compatibility:

- Ensure your smartphone or tablet is compatible with the Dexcom app. Proper pairing with a Wear OS smartwatch relies on having a supported mobile device.

Check Compatibility Requirements:

- For full details, visit Dexcom's official FAQ page to verify requirements and supported devices: [\[Dexcom Compatibility Requirements\]](#).

By understanding these compatibility details, you can confidently use the Dexcom app on your Wear OS smartwatch, ensuring reliable CGM monitoring anytime, anywhere.